


RCCS-4SC


TECHNICAL SPECIFICATIONS-DIMENSIONS RCCS-4SC

| COMPONENT | HEIGHT | WIDTH | DEPTH |
|---------------------|------------------|------------------|-----------------|
| Base | 9.0" (22.9 cm) | 8.0" (20.3 cm) | 8.25" (21.0 cm) |
| Power Supply | 4.5" (11.4 cm) | 9.5" (24.1 cm) | 9.5" (24.1 cm) |
| COMPONENT | WEIGHT | ROTATIONAL SPEED | — |
| Base | 7.0 lbs (3.3 Kg) | 5-27 RPM | — |
| Power Supply | 6.0 lbs (2.83Kg) | — | — |



RCCS-4SC POWER AND ENVIRONMENTAL SPECIFICATIONS

| POWER | VOLTAGE | FREQUENCY |
|----------------------------|-------------------------|-----------|
| Power Supply Input | USA 85-260 VAC @ .5Amps | 50-60 Hz |
| Power Supply Output | 7.5 VDC @ 2Amps | 50-60 Hz |
| ENVIRONMENT | TEMPERATURE | HUMIDITY |
| Base | 50-105 °F (10-40 °C) | 0-100% RH |
| Power Supply | 50-80 °F (10-27 °C) | 0-75% RH |



RCCS-4SC COMPATIBLE CULTURE VESSELS

| CULTURE VESSEL MODEL | SIZE | COMPATIBILITY |
|--------------------------|--------|---------------|
| Disposable HARV * | 10 mL | Yes |
| Disposable HARV * | 50mL | Yes |
| Autoclavable HARV | 1 mL | Yes |
| Autoclavable HARV | 2 mL | Yes |
| Autoclavable HARV | 4 mL | Yes |
| Autoclavable HARV | 10 mL | Yes |
| Autoclavable HARV | 50 mL | No |
| Autoclavable STLV | 55 mL | No |
| Autoclavable STLV | 110 mL | No |
| Autoclavable STLV | 250 mL | No |
| Autoclavable STLV | 500 mL | No |



* Vessels are supplied four to a pack. Four Stem Cell Autoclavable HARVs are provided with the RCCS-4SC System. All other compatible vessel types are ordered separately. Visit www.synthecon.com for more details.

3D Culture Systems
for a 3D World