


HARV

TECHNICAL SPECIFICATIONS-DIMENSIONS HARV (Vessel Chamber- Inner Vessel)

VESSEL	DIAMETER	DEPTH
1mL	0.72" (18.3 mm)	0.15" (3.81 mm)
2mL	.88" (22.23 mm)	0.20" (5.16 mm)
4mL	1.12" (28.32 mm)	0.25" (6.35 mm)
10mL	1.7" (43.18 mm)	0.25" (6.35 mm)
50mL	3.9" (99.06 mm)	0.25" (6.35 mm)



3D Culture Systems
for a 3D World



Primary Business Address
8977 Interchange Road
Houston, Texas
77054

Phone: 713-741-2582
Toll Free: 800-853-0740
Email: sales@synthecon.com
Website: <http://www.synthecon.com>

